

# SAFE.<sup>1</sup> SECURE.<sup>2</sup> EFFECTIVE.<sup>1</sup>

Made possible with the Endo GIA™ reinforced reload with Tri-Staple™ technology<sup>1,2,†</sup>



**You Target Tissue.  
Our Stapler Adjusts Firing Speed.<sup>3,4</sup>**

Adaptive Firing™ technology — built into the iDrive™ Ultra powered stapler — responds to tissue characteristics to optimize staple formation.<sup>3,4</sup>

†As shown in both a clinical and preclinical study.

‡Preclinical results may not correlate with clinical performance in humans.

§Compared to uniform staple heights.

## Study One Overview

A peer-reviewed paper compares Endo GIA™ reinforced reload with Tri-Staple™ technology to Endo GIA™ reload with Tri-Staple™ technology and no buttress.<sup>‡</sup>

## Clinical Challenge

Incomplete staple formation due to variable tissue, displacement of tissue, and staples intersect may lead to factors related to anastomotic failure.<sup>1,‡</sup>

## Technology Solution

Tri-Staple™ technology has been proven to provide superior burst strength.<sup>1,§</sup> When you add a preloaded buttress — the Endo GIA™ reinforced reload with Tri-Staple™ technology — you enhance this performance with improved hemostasis.<sup>2,‡</sup>

## Conclusion

The preloaded buttress on the reinforced reload provides staple formation security when staple lines intersect compared to reloads with no buttress.<sup>2,‡</sup>



## More Consistent Staple Formation


Appropriate formation of B-shape staples is essential for proper tissue apposition.<sup>5</sup>

Complete staple formation rate in colon tissue (thin)<sup>2</sup>



■ Endo GIA™ reinforced reload with Tri-Staple™ technology

Complete staple formation rate in stomach tissue (thick)<sup>2</sup>



■ Endo GIA™ reload with Tri-Staple™ technology

**The rate of complete staple formation was significantly higher with the reinforced reload ( $p < 0.05$ ).**

## High Rate Hypothesis

The authors of this preclinical study<sup>2</sup> attribute the higher rate of complete staple formation to:

- The combination of the Endo GIA™ reinforced reload with Tri-Staple™ technology 60 mm AMT reload with iDrive™ Ultra powered stapling system<sup>2</sup>
- The preloaded material compresses the tissue in the staple line evenly, and prevents tissue displacement when crossing intersecting staple lines<sup>2</sup>

**Medtronic**  
Further, Together



## Secure Staple Formation. Reliable Hemostasis.

The second study showed that the reinforced reload may contribute to reducing — and even eliminating — anastomosis complications when compared to a non-buttressed linear reload.<sup>1</sup>



### 17 PATIENTS

underwent laparoscopic surgery to treat colorectal cancer<sup>1</sup>



### 0 INCIDENTCES

of postoperative bleeding or grade III or higher postoperative complications<sup>1,Ω</sup>

ΩBased on the Clavien–Dindo classification.

## Study Two Overview

### Facilitating a New Technique

According to the second study, when the sigmoid colon, rectosigmoid, or rectum was transected, holding the wings of the material allowed the surgeons to<sup>1</sup>:

- Use the buttressed staple line as a guide to position the trocar of the circular stapler
- Atraumatically grasp the buttress material, instead of the tissue, for circular staple insertion and removal
- Staple the resected rectum in a straight line



### Solving Clinical Challenges

The findings of a careful technique demonstrated how the Endo GIA™ reinforced reload with Tri-Staple™ technology is safe and effective for colorectal resections and anastomoses.<sup>1,2,††</sup>

As the study showed, proper use of the device in gastrointestinal anastomosis can likely reduce complications.<sup>2</sup>

††Results indicate that no major safety issues were reported 30 days post procedure.

1. Naito M, Yamanashi T, Nakamura T, et al. Safety and efficacy of a novel linear staple device with bioabsorbable polyglycolic acid felt in laparoscopic colorectal surgery. *Asian J Endosc Surg*. 2017;10(1):35-39.
2. Naito M, Sato T, Nakamura T, et al. Secure overlap stapling using a linear stapler with bioabsorbable polyglycolic acid felt. *Asian J Endosc Surg*. Feb 22, 2017. [Epub ahead of print]
3. Based on internal test report #R2146-151-0. Powered stapling firing speed DOE analysis and ASA parameters. 2015.
4. Based on internal test report #R2146-173-0. ASA verification testing with slow speed force limit evaluation. 2015.
5. Baker RS, Foote J, Kemmeter P, Brady R, Vroegop T, Serveld M. The science of stapling and leaks. *Obes Surg*. 2004;14:1290–1298.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 04/2018–US170457–[WF#1611073]

**Medtronic**