




ONE COMPLETE SOLUTION.

VersaOne™ Access System






PRODUCT CODE NOMENCLATURE DESIGNED FOR EASY PRODUCT IDENTIFICATION



EXAMPLE

VersaOne™ Blunt Trocar, BPT12STS

BPT = Blunt, 12 = 12 mm,
ST = Standard Length, S = Smooth Cannula



VERSAONE™ ACCESS SYSTEM FULL PORTFOLIO

The VersaOne™ access system streamlines our product codes. This guide provides information on product codes that will change, be discontinued and remain the same.

Description	Code	Size	Length	Units/Box
VersaOne™ Optical Trocars				
VersaOne™ Optical Trocar with Fixation Cannula	ONB5STF	5 mm	100 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB5SHF	5 mm	70 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB5LGF	5 mm	150 mm	6
VersaOne™ Optical Trocar with Dual Fixation Cannula	ONB5STF2C	5 mm	100 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB11STF	11 mm	100 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB11LGF	11 mm	150 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB12STF	12 mm	100 mm	6
VersaOne™ Optical Trocar with Smooth Cannula	ONB12STS	12 mm	100 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB12SHF	12 mm	70 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB12LGF	12 mm	150 mm	6
VersaOne™ Optical Trocar with Fixation Cannula	ONB15STF	15 mm	100 mm	6
VersaOne™ Blunt Trocar				
VersaOne™ Blunt Trocar with Threaded Anchor	BPT12STS	12 mm	100 mm	3



Description	Code	Size	Length	Units/Box
VersaOne™ Bladed Trocars				
VersaOne™ Bladed Trocar with Fixation Cannula	B5STF	5 mm	100 mm	6
VersaOne™ Bladed Trocar with Smooth Cannula	B5STS	5 mm	100 mm	6
VersaOne™ Bladed Trocar with Fixation Cannula	B5SHF	5 mm	70 mm	6
VersaOne™ Bladed Trocar with Smooth Cannula	B5SHS	5 mm	70 mm	6
VersaOne™ Bladed Trocar with Fixation Cannula	B11STF	11 mm	100 mm	6
VersaOne™ Bladed Trocar with Smooth Cannula	B11STS	11 mm	100 mm	6
VersaOne™ Bladed Trocar with Fixation Cannula	B12STF	12 mm	100 mm	6
VersaOne™ Bladed Trocar with Smooth Cannula	B12STS	12 mm	100 mm	6
VersaOne™ Bladed Trocar with Fixation Cannula	B12LGF	12 mm	150 mm	6
VersaOne™ Bladed Trocar with Smooth Cannula	B12LGS	12 mm	150 mm	6
VersaOne™ Bladeless Trocars				
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5STF	5 mm	100 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5SHF	5 mm	70 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5LGF	5 mm	150 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB8STF	8 mm	100 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB11STF	11 mm	100 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB12STF	12 mm	100 mm	6
VersaOne™ Bladeless Trocar with Smooth Cannula	NONB12STS	12 mm	100 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB12SHF	12 mm	70 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB12LGF	12 mm	150 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB15STF	15 mm	100 mm	6
VersaOne™ Bladeless Trocar with Fixation Cannula	NONB15LGF	15 mm	150 mm	6
VersaOne™ Universal Cannula				
VersaOne™ Universal Fixation Cannula	UNVCA5STF	5 mm	100 mm	6
VersaOne™ Universal Fixation Cannula	UNVCA5SHF	5 mm	70 mm	6
VersaOne™ Universal Fixation Cannula	UNVCA5LGF	5 mm	150 mm	6
VersaOne™ Universal Fixation Cannula	UNVCA11STF	11 mm	100 mm	6
VersaOne™ Universal Fixation Cannula	UNVCA12STF	12 mm	100 mm	6
VersaOne™ Universal Smooth Cannula	UNVCA12STS	12 mm	100 mm	6
VersaOne™ Universal Fixation Cannula	UNVCA12LGF	12 mm	150 mm	6



LEGACY PRODUCT CODE UPDATE

The VersaOne™ access system streamlines our product codes. The following pages identify the changes being made to codes, packaging and trocar diameters.

Legacy Description	Legacy Code	New Description	New Code	Length	Units/Box
Versaport™ Bladeless Optical Trocars		VersaOne™ Optical Trocars			
Versaport™ Bladeless Optical 5 mm Standard Fixation	ONB5STF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	100 mm	6
Versaport™ Bladeless Optical 5 mm Short Fixation	ONB5SHF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	70 mm	6
Versaport™ Bladeless Optical 5 mm Long Fixation	ONB5LGF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	150 mm	6
Versaport™ Bladeless Optical 5 mm Standard Fixation, Dual Pack	ONB5STF2C	VersaOne™ Optical Trocar with Fixation Cannula	No Change	100 mm	6
Versaport™ Bladeless Optical 11 mm Standard Fixation	ONB11STF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	100 mm	6
Versaport™ Bladeless Optical 11 mm Long Fixation	ONB11LGF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	150 mm	6
Versaport™ Bladeless Optical 12 mm Standard Fixation	ONB12STF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	100 mm	6
Versaport™ Bladeless Optical 12 mm Standard Smooth	ONB12STS	VersaOne™ Optical Trocar with Smooth Cannula	No Change	100 mm	6
Versaport™ Bladeless Optical 12 mm Short Fixation	ONB12SHF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	70 mm	6
Versaport™ Bladeless Optical 12 mm Long Fixation	ONB12LGF	VersaOne™ Optical Trocar with Fixation Cannula	No Change	150 mm	6
Versaport™ Bladeless Optical 5 mm Standard Fixation, Cannula Only	ONBFCA5ST	VersaOne™ Universal Fixation Cannula	UNVCA5STF	100 mm	6
Versaport™ Bladeless Optical 5 mm Short Fixation, Cannula Only	ONBFCA5SH	VersaOne™ Universal Fixation Cannula	UNVCA5SHF	70 mm	6
Versaport™ Bladeless Optical 5 mm Long Fixation, Cannula Only	ONBFCA5LG	VersaOne™ Universal Fixation Cannula	UNVCA5LGF	150 mm	6
Versaport™ Bladeless Optical 11 mm Standard Fixation, Cannula Only	ONBFCA11ST	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ Bladeless Optical 12 mm Standard Fixation, Cannula Only	ONBFCA12ST	VersaOne™ Universal Fixation Cannula	UNVCA12STF	100 mm	6
-----	-----	VersaOne™ Universal Smooth Cannula	UNVCA12STS	100 mm	6
Versaport™ Bladeless Optical 12 mm Long Fixation, Cannula Only	ONBFCA12LG	VersaOne™ Universal Fixation Cannula	UNVCA12LGF	150 mm	6
Bluntport™ Plus Trocars*		VersaOne™ Blunt Trocars			
Bluntport™ Plus 12 mm with FoamGrip™ Standard Smooth	176626P	VersaOne™ Blunt Trocar with Threaded Anchor	BPT12STS	100 mm	3
Bluntport™ Plus 12 mm with BluntGrip™ Standard Smooth	179075P	VersaOne™ Blunt Trocar with Threaded Anchor	BPT12STS	100 mm	3



Legacy Description	Legacy Code	New Description	New Code	Length	Units/Box
Versaport™ Bladed V2 Trocars*		VersaOne™ Bladed Trocars			
Versaport™ V2 5 mm Standard Fixation	179094F	VersaOne™ Bladed Trocar with Fixation Cannula	B5STF	100 mm	6
Versaport™ V2 5 mm Standard Smooth	179094	VersaOne™ Bladed Trocar with Smooth Cannula	B5STS	100 mm	6
Versaport™ V2 5 mm Short Fixation	179093F	VersaOne™ Bladed Trocar with Fixation Cannula	B5SHF	70 mm	6
Versaport™ V2 5 mm Short Smooth	179093	VersaOne™ Bladed Trocar with Smooth Cannula	B5SHS	70 mm	6
Versaport™ V2 11 mm Standard Fixation	179095PF	VersaOne™ Bladed Trocar with Fixation Cannula	B11STF	100 mm	6
Versaport™ V2 11 mm Standard Smooth	179095P	VersaOne™ Bladed Trocar with Smooth Cannula	B11STS	100 mm	6
Versaport™ V2 12 mm Standard Fixation	179096PF	VersaOne™ Bladed Trocar with Fixation Cannula	B12STF	100 mm	6
Versaport™ V2 12 mm Standard Smooth	179096P	VersaOne™ Bladed Trocar with Smooth Cannula	B12STS	100 mm	6
Versaport™ V2 12 mm Long Fixation	179097PF	VersaOne™ Bladed Trocar with Fixation Cannula	B12LGF	150 mm	6
Versaport™ V2 12 mm Long Smooth	179097P	VersaOne™ Bladed Trocar with Smooth Cannula	B12LGS	150 mm	6
Versaport™ V2 5 mm Standard Fixation, Cannula Only	177092F	VersaOne™ Universal Fixation Cannula	UNVCA5STF	100 mm	6
Versaport™ V2 5 mm Standard Smooth, Cannula Only	177092	VersaOne™ Universal Fixation Cannula	UNVCA5STF	100 mm	6
Versaport™ V2 5 mm Standard Smooth, Cannula Only	177990	VersaOne™ Universal Fixation Cannula	UNVCA5STF	100 mm	6
Versaport™ V2 5 mm Short Fixation, Cannula Only	177093F	VersaOne™ Universal Fixation Cannula	UNVCA5SHF	70 mm	6
Versaport™ V2 5 mm Short Smooth, Cannula Only	177093	VersaOne™ Universal Fixation Cannula	UNVCA5SHF	70 mm	6
Versaport™ V2 5 mm Short Smooth, Cannula Only	177991	VersaOne™ Universal Fixation Cannula	UNVCA5SHF	70 mm	6
Versaport™ V2 10 mm Standard Smooth, Cannula Only	177096P	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ V2 10 mm Standard Smooth, Cannula Only	177994P	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ V2 11 mm Standard Fixation, Cannula Only	177090PF	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ V2 11 mm Standard Smooth, Cannula Only	177090P	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ V2 11 mm Standard Smooth, Cannula Only	177992P	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ V2 12 mm Standard Fixation, Cannula Only	177091PF	VersaOne™ Universal Fixation Cannula	UNVCA12STF	100 mm	6
Versaport™ V2 12 mm Standard Smooth, Cannula Only	177995P	VersaOne™ Universal Smooth Cannula	UNVCA12STS	100 mm	6



*Previously a box quantity of 3

LEGACY PRODUCT CODE UPDATE

The VersaOne™ access system streamlines our product codes. This guide provides information on product codes that will change, be discontinued and remain the same.

Legacy Description	Legacy Code	New Description	New Code	Length	Units/Box
Versaport™ RPF Bladed Trocars		VersaOne™ Bladed Trocars			
Versaport™ Plus RPF Bladed 8 mm Standard Smooth	179074P	No Replacement / To Be Discontinued			
Versaport™ Plus RPF Bladed 10 mm Standard Smooth	179076P	No Replacement / To Be Discontinued			
Versaport™ Plus RPF Bladed 10 mm Long Smooth	179077P	No Replacement / To Be Discontinued			
Versaport™ Plus RPF 15 mm Standard Smooth	179078P	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB15STF	100 mm	6
		VersaOne™ Optical Trocar with Fixation Cannula	ONB15STF		
Versaport™ Plus RPF 15 mm Standard Fixation	179078PF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB15STF	100 mm	6
		VersaOne™ Optical Trocar with Fixation Cannula	ONB15STF		
Dexide™ Trocar		VersaOne™ Bladed Trocars			
Dexide™ 5 mm Threaded Trocar	4016	VersaOne™ Bladed Trocar with Fixation Cannula	B5SHF	70 mm	6
Dexide™ 5 mm Threaded Trocar w/2 Cannulas	4216	VersaOne™ Bladed Trocar with Fixation Cannula	B5SHF	70 mm	6
		VersaOne™ Universal Fixation Cannula	UNVCA5SHF	70 mm	6



Legacy Description	Legacy Code	New Description	New Code	Length	Units/Box
Versaport™ Bladeless Trocars		VersaOne™ Bladeless Trocars			
Versaport™ Bladeless Low Profile 2/3 mm Short Fixation	NB3SHFLP	No Change	No Change	70 mm	6
Versaport™ Bladeless Low Profile 5 mm Standard Fixation	NB5STFLP	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5STF	100 mm	6
Versaport™ Bladeless Low Profile 5 mm Short Fixation	NB5SHFLP	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5SHF	70 mm	6
Versaport™ Bladeless 5 mm Standard Fixation	NB5STF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5STF	100 mm	6
Versaport™ Bladeless 5 mm Short Fixation	NB5SHF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5SHF	70 mm	6
Versaport™ Plus Bladeless 5 mm Long Fixation	NB5LGF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB5LGF	150 mm	6
Versaport™ Bladeless 8 mm Standard Fixation	NB8STF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB8STF	100 mm	6
Versaport™ Plus Bladeless 11 mm Standard Fixation	NB11STF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB11STF	100 mm	6
Versaport™ Plus Bladeless 11 mm Standard Smooth	NB11STS	No Replacement / To Be Discontinued			
Versaport™ Plus Bladeless 12 mm Standard Fixation	NB12STF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB12STF	100 mm	6
Versaport™ Plus Bladeless 12 mm Standard Smooth	NB12STS	VersaOne™ Bladeless Trocar with Smooth Cannula	NONB12STS	100 mm	6
Versaport™ Plus Bladeless 12 mm Short Fixation	NB12SHF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB12SHF	70 mm	6
Versaport™ Plus Bladeless 12 mm Long Fixation	NB12LGF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB12LGF	150 mm	6
Versaport™ Plus Bladeless 15 mm Standard Fixation	NB15STF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB15STF	100 mm	6
Versaport™ Plus Bladeless 15 mm Long Fixation	NB15LGF	VersaOne™ Bladeless Trocar with Fixation Cannula	NONB15LGF	150 mm	6
Versaport™ Plus 5 mm Standard Fixation, Cannula Only	NBFCA5ST	VersaOne™ Universal Fixation Cannula	UNVCA5STF	100 mm	6
Versaport™ Plus 5 mm Short Fixation, Cannula Only	NBFCA5SH	VersaOne™ Universal Fixation Cannula	UNVCA5SHF	70 mm	6
Versaport™ Plus 11 mm Standard Fixation, Cannula Only	NBFCA11ST	VersaOne™ Universal Fixation Cannula	UNVCA11STF	100 mm	6
Versaport™ Plus 12 mm Standard Fixation, Cannula Only	NBFCA12ST	VersaOne™ Universal Fixation Cannula	UNVCA12STF	100 mm	6



VERSAONE™ ACCESS SYSTEM PRODUCT CODE MEASUREMENTS

VersaOne™ Access System Product Code Measurements

VersaOne™ Access System Products	Inner Diameter	Outer Diameter
VersaOne™ Access System 5 mm	5.92 mm	8.46 mm
VersaOne™ Access System 8 mm	8.64 mm	11.18 mm
VersaOne™ Access System 11 mm	11.47 mm	13.97 mm
VersaOne™ Access System 12 mm	12.95 mm	15.52 mm
VersaOne™ Access System 15 mm	15.52 mm	18.02 mm

Legacy Product Code Measurements

Versaport™ Optical	Inner Diameter	Outer Diameter
5 mm	5.92 mm	8.46 mm
11 mm	11.47 mm	13.97 mm
12 mm	12.95 mm	15.52 mm

Versaport™ V2 Bladed	Inner Diameter	Outer Diameter
5 mm	5.56 mm	8.00 mm
11 mm	11.33 mm	13.97 mm
12 mm	12.45 mm	15.14 mm

Versaport™ Plus Bladeless	Inner Diameter	Outer Diameter
5 mm	5.92 mm	8.20 mm
8 mm	8.38 mm	11.07 mm
11 mm	11.33 mm	13.97 mm
12 mm	12.95 mm	15.52 mm
15 mm	15.47 mm	18.14 mm

Bluntport™ Plus	Inner Diameter	Outer Diameter
12 mm	12.45 mm	14.53 mm





Contact your Medtronic representative
for more information at
800-722-8772 or visit our website
medtronic.com/covidien

555 Long Wharf Drive
New Haven, CT 06511

Phone: 800.722.8772
Fax: 888.636.1002

[medtronic.com/covidien](https://www.medtronic.com/covidien)