

YOUR GUIDE  
TO OUR  
**CARDIOVASCULAR  
WOUND CLOSURE  
SOLUTIONS.**



**Cardiovascular wound closure  
cross reference guide**

**Medtronic**  
Further, Together

# NEW CODE GUIDE

## SUTURES

### Flexon™ Suture

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
8886258643	8886258643-2	FLXN 3-0 24 CL CV-23 SC-6 DA	Size 3-0	CV-23	SC-6	1 Clear	24	60	-	Regular Attach	1	12	Taper Point	Conventional Cutting	Cardiopoint™ Needle
8886259143	8886259143-2	FLXN 3-0 24 CL CV-23 SC DA	Size 3-0	CV-23	SC-6	1 Clear	24	60	-	Regular Attach	1	12	Taper Point	Conventional Cutting	Cardiopoint™ Needle
8886259243	8886259243-2	FLXN 3-0 24 CL CV-23 SC-2 DA	Size 3-0	CV-23	SC-2	1 Clear	24	60	-	Regular Attach	1	12	Taper Point	Conventional Cutting	Cardiopoint™ Needle
8886259363	8886259363-2	FLXN 0 24 CL CV-24 SC-2 DA	Size 0	CV-24	SC-2	5 Blue, 5 White	24	60	-	Regular Attach	1	12	Taper Point	Conventional Cutting	Cardiopoint™ Needle

### Ti-Cron™ Suture

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
8886263656	8886263656-2	TICRON 2-0 10X36 BL-WH CV-15 DA	Size 2-0	CV-15	CV-15	5 Blue, 5 White	36	90	-	Regular Attach	5	6	Taper Point	Taper Point	Cardiopoint™ Needle
8886310421	8886310421-2	TICRON 5-0 30 BL CV-11	Size 5-0	CV-11	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886326742	8886326742-2	TICRON 3-0 18 BL CV-24	Size 3-0	CV-24	-	1 Blue	18	45	-	D-Tach	5	24	Taper Point	-	Cardiopoint™ Needle

### Maxon™ Suture

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
8886603641	8886603641-2	MAXON 3-0 36 GRN CV-24 DA	Size 3-0	CV-24	CV-24	1 Green	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886605731	8886605731-2	MAXON CV 4-0 30 GRN CV-15	Size 4-0	CV-15	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886612111	8886612111-2	MAXON CV 6-0 18 GRN CV-11 DA	Size 6-0	CV-11	CV-11	1 Green	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886612411	8886612411-2	MAXON CV 6-0 30 GRN CV-11 DA	Size 6-0	CV-11	CV-11	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886612421	8886612421-2	MAXON CV 5-0 30 GRN CV-11 DA	Size 5-0	CV-11	CV-11	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886612711	8886612711-2	MAXON CV 6-0 30 GRN CV-1 DA	Size 6-0	CV-1	CV-1	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886612821	8886612821-2	MAXON CV 5-0 36 GRN CV-11 DA	Size 5-0	CV-11	CV-11	1 Green	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886613431	8886613431-2	MAXON CV 4-0 36 GRN CV-23 DA	Size 4-0	CV-23	CV-23	1 Green	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886613911	8886613911-2	MAXON CV 6-0 24 GR CV-1 DA	Size 6-0	CV-1	CV-1	1 Green	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886614711	8886614711-2	MAXON CV 6-0 24 GRN CV-11 DA	Size 6-0	CV-11	CV-11	1 Green	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886614721	8886614721-2	MAXON CV 5-0 24 GRN CV-11 DA	Size 5-0	CV-11	CV-11	1 Green	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886614811	8886614811-2	MAXON CV 6-0 30 GRN CV-11	Size 6-0	CV-11	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886614821	8886614821-2	MAXON CV 5-0 30 GRN CV-301 CV-11	Size 5-0	CV-11	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886614921	8886614921-2	MAXON CV 5-0 30 GRN CV-23 DA	Size 5-0	CV-23	CV-23	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886620721	8886620721-2	MAXON CV 5-0 30 GRN CV-23	Size 5-0	CV-23	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886620731	8886620731-2	MAXON CV 4-0 30 GRN CV-23	Size 4-0	CV-23	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886620741	8886620741-2	MAXON CV 3-0 30 GRN CV-23	Size 3-0	CV-23	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886621721	8886621721-2	MAXON CV 5-0 30 GRN CV-24	Size 5-0	CV-24	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886621731	8886621731-2	MAXON CV 4-0 30 GRN CV-24	Size 4-0	CV-24	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886621741	8886621741-2	MAXON 3-0 30 GRN CV-24	Size 3-0	CV-24	-	1 Green	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886621943	8886621943-2	MAXON 3-0 18 GRN CV-24	Size 3-0	CV-24	-	1 Green	18	45	-	D-Tach	5	12	Taper Point	-	Cardiopoint™ Needle
8886622131	8886622131-2	MAXON CV 4-0 36 GRN CV-24 DA	Size 4-0	CV-24	CV-24	1 Green	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886646221	8886646221-2	MAXON CV 5-0 30 GRN CV-15 DA	Size 5-0	CV-15	CV-15	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886646231	8886646231-2	MAXON CV 4-0 30 GRN CV-15 DA	Size 4-0	CV-15	CV-15	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
XX2121	XX2121-2	MAXON CV 6-0 18 CLEAR CV-1	Size 6-0	CV-1	-	1 Clear	18	45	-	Regular Attach	1	12	Taper Point	-	Cardiopoint™ Needle
8886612701	8886612701-2	MAXON CV 7-0 30 GRN CV-1 DA	Size 7-0	CV-1	CV-1	1 Green	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886613901	8886613901-2	MAXON CV 7-0 24 GRN CV-1 DA	Size 7-0	CV-1	CV-1	1 Green	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886614501	8886614501-2	MAXON CV 7-0 18 GRN CV-1 DA	Size 7-0	CV-1	CV-1	1 Green	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886616523	8886616523-2	MAXON CV 5-0 4X24 GRN CV-11 DA	Size 5-0	CV-11	CV-11	1 Green	24	60	-	Regular Attach	4	12	Taper Point	Taper Point	Cardiopoint™ Needle
8886621933	8886621933-2	MAXON CV 4-0 5X18 GRN CV-24 DT	Size 4-0	CV-24	-	1 Green	18	45	-	D-Tach	5	12	Taper Point	-	Cardiopoint™ Needle

# NEW CODE GUIDE

## SUTURES

### Sofsilkm Suture

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
GS451	GS451-2	SOFILK 4-0 30 BK CV-24	Size 4-0	CV-24	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
GS452	GS452-2	SOFILK 3-0 30 BK CV-24	Size 3-0	CV-24	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
GS453	GS453-2	SOFILK 2-0 30 BK CV-24	Size 2-0	CV-24	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
GS46M	GS46M-2	SOFILK 4-0 5x18 BK CV-23 DT	Size 4-0	CV-23	-	1 Black	18	45	-	D-Tach	5	24	Taper Point	-	Cardioint™ Needle
GS47M	GS47M-2	SOFILK 3-0 5x18 BK CV-23 DT	Size 3-0	CV-23	-	1 Black	18	45	-	D-Tach	5	24	Taper Point	-	Cardioint™ Needle
GSJ36M	GSJ36M-2	SOFSLK 4-0 5X18 BK CV-24 DT	Size 4-0	CV-24	-	1 Black	18	45	-	D-Tach	5	24	Taper Point	-	Cardioint™ Needle
GSJ37M	GSJ37M-2	SOFILK 3-0 5X18 BK CV-24 DT	Size 3-0	CV-24	-	1 Black	18	45	-	D-Tach	5	24	Taper Point	-	Cardioint™ Needle
GSJ46M	GSJ46M-2	SOFILK 4-0 5x18 BK CV-23 DT	Size 4-0	CV-23	-	1 Black	18	45	-	D-Tach	5	24	Taper Point	-	Cardioint™ Needle
GSJ47M	GSJ47M-2	SOFILK 3-0 5x18 BK CV-23 DT	Size 3-0	CV-23	-	1 Black	18	45	-	D-Tach	5	24	Taper Point	-	Cardioint™ Needle
VS581	VS581-2	SOFILK 4-0 24 BK CV-15 DA	Size 4-0	CV-15	CV-15	1 Black	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VS706M	VS706M-2	SOFILK 6-0 5x24 BK CV-11	Size 6-0	CV-11	CV-11	1 Black	24	60	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VS709	VS709-2	SOFILK 4-0 18 BK CV-11 DA	Size 4-0	CV-11	CV-11	1 Black	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VS802	VS802-2	SOFILK 6-0 18 BK CV-1	Size 6-0	CV-1	-	1 Black	18	45	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS806	VS806-2	SOFILK 6-0 18 BK CV-11 DA	Size 6-0	CV-11	CV-11	1 Black	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VS809	VS809-2	SOFILK 7-0 18 BK CV-1 DA	Size 7-0	CV-1	CV-1	1 Black	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VS823	VS823-2	SOFILK 2-0 30 BL CV-15M	Size 2-0	CV-15	CV-15	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS870	VS870-2	SOFILK 5-0 30 BK CV-23	Size 5-0	CV-23	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS871	VS871-2	SOFILK 4-0 30 BK CV-23	Size 4-0	CV-23	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS872	VS872-2	SOFILK 3-0 30 BK CV-23	Size 3-0	CV-23	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS873	VS873-2	SOFILK 2-0 30 BK CV-23	Size 2-0	CV-23	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS880	VS880-2	SOFILK 5-0 30 BLK CV-15	Size 5-0	CV-15	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS881	VS881-2	SOFILK 4-0 30 BK CV-15	Size 4-0	CV-15	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS882	VS882-2	SOFILK 3-0 30 BK CV-15	Size 3-0	CV-15	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS889	VS889-2	SOFILK 6-0 30 BK CV-11	Size 6-0	CV-11	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS890	VS890-2	SOFILK 5-0 30 BK CV-11	Size 5-0	CV-11	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS891	VS891-2	SOFILK 4-0 30 BK CV-11	Size 4-0	CV-11	-	1 Black	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VS89G	VS89G-2	SOFILK 7-0 18 BL CV-1M DA	Size 7-0	CV-1	CV-1	1 Black	18	45	-	Regular Attach	12	12	Taper Point	Taper Point	Cardioint™ Needle
VS823	VS823-2	SOFILK 2-0 30 BK CV-15	Size 2-0	CV-15	-	1 Black	30	75	-	Reg. Attach			Taper Point	-	Cardioint™ Needle

### Surgipro™ II Suture

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
VP331X	VP331X-2	SRGPROII 4-0 36 BL CV-24 DA	Size 4-0	CV-24	CV-24	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP363X	VP363X-2	SRGPROII 6-0 36 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP365X	VP365X-2	SRGPROII 6-0 36 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP555X	VP555X-2	SRGPROII 6-0 36 BLCV-23 DA	Size 6-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP556MX	VP556MX-2	SRGPROII 5-0 5X36 BL CV-23 DA	Size 5-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP556X	VP556X-2	SRGPROII 5-0 36 BL CV-23 DA	Size 5-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP557MX	VP557MX-2	SRGPROII 4-0 5X36 BL CV-23 DA	Size 4-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP557X	VP557X-2	SRGPROII 4-0 36 BL CV-23 DA	Size 4-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP558X	VP558X-2	SRGPROII 3-0 36 BL CV-23 DA	Size 3-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP566X	VP566X-2	SRGPROII 5-0 24 BL CV-23 DA	Size 5-0	CV-23	CV-23	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP570X	VP570X-2	SRGPROII 5-0 30 BL CV-15 DA	Size 5-0	CV-15	CV-15	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP580X	VP580X-2	SRGPROII 5-0 36 BL CV-15 DA	Size 5-0	CV-15	CV-15	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP581X	VP581X-2	SRGPROII 4-0 36 BL CV-15 DA	Size 4-0	CV-15	CV-15	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP582X	VP582X-2	SRGPROII 4-0 36 BL CV-15 DA	Size 4-0	CV-15	CV-15	1 Blue	36	90	3 x 7 mm Firm	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP583X	VP583X-2	SRGPROII 4-0 48 BL CV-15 DA	Size 4-0	CV-15	CV-15	1 Blue	48	120	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP584X	VP584X-2	SRGPROII 4-0 24 BL CV-15 DA	Size 4-0	CV-15	CV-15	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP585X	VP585X-2	SRGPRO II 3-0 36 BL CV-15 DA	Size 3-0	CV-15	CV-15	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP701MX	VP701MX-2	SRGPROII 7-0 5X18 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	18	45	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP701X	VP701X-2	SRGPROII 7-0 18 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP702MX	VP702MX-2	SRGPROII 7-0 5X24 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle

# NEW CODE GUIDE

## SUTURES

### Surgipro™ II Suture cont.

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
VP702X	VP702X-2	SRGPROII 7-0 24 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP703MX	VP703MX-2	SRGPROII 6-0 5X24 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP703X	VP703X-2	SRGPROII 6-0 24 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP705MX	VP705MX-2	SRGPROII 5-0 5X30 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP705X	VP705X-2	SRGPROII 5-0 30 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP706MX	VP706MX-2	SRGPROII 6-0 5X30 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP706X	VP706X-2	SRGPROII 6-0 30 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP709MX	VP709MX-2	SRGPROII 6-0 5X30 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP709X	VP709X-2	SRGPROII 6-0 30 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP713MX	VP713MX-2	SRGPROII 7-0 5X30 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP713X	VP713X-2	SRGPROII 7-0 30 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP718MX	VP718MX-2	SRGPROII 6-0 5X18 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP718X	VP718X-2	SRGPROII 6-0 18 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP719MX	VP719MX-2	SRGPROII 5-0 5X18 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP719X	VP719X-2	SRGPROII 5-0 18 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP720MX	VP720MX-2	SRGPROII 5-0 5X36 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	36	90	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP720X	VP720X-2	SRGPROII 5-0 36 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP725MX	VP725MX-2	SRGPROII 5-0 5X24 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP725X	VP725X-2	SRGPROII 5-0 24 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP726MX	VP726MX-2	SRGPROII 6-0 5X24 BL CV11 DA	Size 6-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	5	12	Taper Point	Taper Point	Cardioint™ Needle
VP726X	VP726X-2	SRGPRO II 6-0 24 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP727X	VP727X-2	SRGPROII 7-0 24 BL CV-11 DA	Size 7-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP746MX	VP746MX-2	SRGPROII 6-0 2x18 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	18	45	-	Regular Attach	2	12	Taper Point	Taper Point	Cardioint™ Needle
VP746X	VP746X-2	SRGPROII 6-0 18 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP757X	VP757X-2	SRGPROII 4-0 18 BL CV-23 DA	Size 4-0	CV-23	CV-23	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP800X	VP800X-2	SRGPROII 7-0 18 BL CV-11 DA	Size 7-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP813X	VP813X-2	SRGPROII 4-0 30 BL CV-11 DA	Size 4-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP870X	VP870X-2	SRGPROII 5-0 30 BL CV-23	Size 5-0	CV-23	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP871X	VP871X-2	SRGPROII 4-0 30 BL CV-23	Size 4-0	CV-23	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP872X	VP872X-2	SRGPROII 3-0 30 BL CV-23	Size 3-0	CV-23	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP881X	VP881X-2	SRGPROII 4-0 30 BL CV-15	Size 4-0	CV-15	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP889X	VP889X-2	SRGPRO II 6-0 30 BL CV-11	Size 6-0	CV-11	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP890X	VP890X-2	SRGPROII 5-0 30 BL CV-11	Size 5-0	CV-11	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP891	VP891-2	SRGPRO 4-0 30 BL CV-11	Size 4-0	CV-11	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardioint™ Needle
VP897M	VP897M-2	SRGPRO 5-0 2x30 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	2	12	Taper Point	Taper Point	Cardioint™ Needle
VPP557X	VPP557X-2	SRGPROII 4-0 36 BL CV23 DA 3X7F	Size 4-0	CV-23	CV-23	1 Blue	36	90	3 x 7 mm Firm	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP559X	VP559X-2	SRGPROII 2-0 BL 36 CV-23M DA	Size 2-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP660X	VP660X-2	SRGPROII 5-0 36 BL CV-23M DA	Size 5-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP78X	VP78X-2	SRGPROII 6-0 18 BL CV-11M DA	Size 6-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle
VP812X	VP812X-2	SRGPROII 7-0 30 BL CV-11M DA	Size 7-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardioint™ Needle

# NEW CODE GUIDE

## SUTURES

### Vascufil™ Suture

Old Existing Codes	New Code	Description	USP Size	Needle 1	Needle 2	Strand Color	Length inches	Length cms	Pledgets	Needle Attachment	Strand Per Box	Box Quantity	Needle Type 1	Needle Type 2	Comments
8886470105V	8886470105V-2	VSCUFL 7-0 4X30 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	4	12	Taper Point	Taper Point	Cardiopoint™ Needle
8886470205V	8886470205V-2	VSCUFL 7-0 4X24 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	4	12	Taper Point	Taper Point	Cardiopoint™ Needle
8886470301V	8886470301V-2	VSCUFL 7-0 18 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470311V	8886470311V-2	VSCUFL 6-0 18 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470401V	8886470401V-2	VSCUFL 7-0 30 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470411V	8886470411V-2	VSCUFL 6-0 30 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470501V	8886470501V-2	VSCUFL 7-0 24 BL CV-1 DA	Size 7-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470511V	8886470511V-2	VSCUFL 6-0 24 BL CV-1 DA	Size 6-0	CV-1	CV-1	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470711V	8886470711V-2	VSCUFL 6-0 18 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	18	45	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470911V	8886470911V-2	VSCUFL 6-0 24 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886470921V	8886470921V-2	VSCUFL 5-0 24 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	24	60	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471011V	8886471011V-2	VSCUFL 6-0 30 BL CV-11 DA	Size 6-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471021V	8886471021V-2	VSCUFL 5-0 30 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471121V	8886471121V-2	VSCUFL 5-0 36 BL CV-11 DA	Size 5-0	CV-11	CV-11	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471621V	8886471621V-2	VSCUFL 5-0 36 BL CV-15 DA	Size 5-0	CV-15	CV-15	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471631V	8886471631V-2	VSCUFL 4-0 36 BL CV-15 DA	Size 4-0	CV-15	CV-15	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471921V	8886471921V-2	VSCUFL 5-0 36 BL CV-23 DA	Size 5-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886471931V	8886471931V-2	VSCUFL 4-0 36 BL CV-23 DA	Size 4-0	CV-23	CV-23	1 Blue	36	90	-	Regular Attach	1	36	Taper Point	Taper Point	Cardiopoint™ Needle
8886472221V	8886472221V-2	VSCUFL 5-0 30 BL CV-11	Size 5-0	CV-11	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886473921V	8886473921V-2	VSCUFL 5-0 30 BL CV-15	Size 5-0	CV-15	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle
8886473931V	8886473931V-2	VSCUFL 4-0 30 BL CV-15	Size 4-0	CV-15	-	1 Blue	30	75	-	Regular Attach	1	36	Taper Point	-	Cardiopoint™ Needle

## HOW TO ORDER

Contact your Medtronic representative for more information.

**800-722-8772**

Visit us at:

**[medtronic.com/covidien](https://www.medtronic.com/covidien)**

555 Long Wharf Drive  
New Haven, CT 06511  
Phone: 1-800-722-8772  
Fax: 1-888-636-1002