





Product Name	Description	Order Code	Box Quantity
ReliaTack™ Articulating Reloadable Fixation Device  <p>with Deep Purchase Tacks One handle with 3 x 8-Tack Deep Purchase Reloads and 1 x 5-Tack Deep Purchase Reloads</p>		RELTACK4XDPT	6
	Eight-tack deep purchase reload	RELTACK8RDPT	6
	Five-tack deep purchase reload	RELTACK5RDPT	6
ReliaTack™ Articulating Reloadable Fixation Device  <p>with Standard Purchase Tacks One handle with 3 x 10-Tack Standard Purchase Reloads</p>		RELTACK3X10	6
	Ten-tack standard purchase reload	RELTACK10R	6
	Five-tack standard purchase reload	RELTACK5R	6
AbsorbaTack™ Absorbable Fixation Device  <p>5 mm Single Use Absorbable Synthetic Polyester Copolymer 30 Tacks Long with Flex Cable Shaft for Laparoscopic Hernia Repair</p> <p>Tissue Purchase Depth 4.1 mm Proximal Wings 5.1 mm</p>		ABSTACK30X	6
	5 mm Single Use Absorbable Synthetic PGLA 30 Tacks Long for Laparoscopic Hernia Repair	ABSTACK30	6
	5 mm Single Use Absorbable Synthetic PGLA 20 Tacks Long for Open Inguinal and Open Ventral Hernia Repair	ABSTACK20S	6
	5 mm Single Use Absorbable Synthetic Polyester Copolymer 15 Tacks Short for Laparoscopic Inguinal Hernia Repair	ABSTACK15	6
ProTack™ Titanium Fixation Device  <p>5 mm Single Use Instrument with 30 Titanium Helical Fasteners with Low Profile Titanium Mesh Interface</p> <p>Thread Width 3.96 mm</p>		174006	6

For more information, call **800-772-8772**
or visit us at [medtronic.com/covidien](https://www.medtronic.com/covidien)